A WIDE RANGE TO MEET ALL YOUR NEEDS









A WIDE RANGE TO MEET ALL YOUR NEEDS













FUNNEL WITH HINGED LID FOR RESIDUAL SOLVENTS

It is installed on a can or a barrel of residual solvent. It is practical and you can equip it with an activated charcoal cartridge if the storage is in a non-ventilated room. The activated charcoal cartridge prevents the creation of an overpressure in the can if it is exposed to heat.



DESCRIPTION	REFERENCE WITH CARTRIDGE*	REFERENCE WITH- OUT CARTRIDGE
Smart Waste Caps funnel with hinged lid for can S55	SW55-FSC	SW55-FSB
Smart Waste Caps funnel with hinged lid for can S60	SW60-FSC	SW60-FSB
Smart Waste Caps funnel with hinged lid for can S90	SW90-FSC	SW90-FSB
Smart Waste Caps funnel with hinged lid for space-saving can St	50 SWB-FSC	SWB-FSB
Smart Waste Caps funnel with hinged lid for metallic barrel	SWBA-FSC	SWBA-FSB

^{*}Activated charcoal cartridge sold separately (see page 24)


Funnel with hinged lid for residual solvents

This funnel has a double interest: the direct discharge through the funnel of all your waste, but also the collection of effluent from your HPLC systems with the integration of inlet fittings and / or threaded barbed fittings on the funnel.

DESCRIPTION	THREAD	REFERENCE
	S55	SW55-2-1L-FSC
Funnel with hinged lid with 1	S60	SW60-2-1L-FSC
activated charcoal cartridge, 2 inlets 1/8" (3.2mm) and 1	S90	SW90-2-1L-FSC
threaded barbed fitting (6 to 9mm)	Metallic barrel	SWBA-2-1L-FSC
	Space-saving can	SWB-2-1L-FSC
	S55	SW55-2-2L-FSC
Funnel with hinged lid with 1	S60	SW60-2-2L-FSC
activated charcoal cartridge, 2 inlets 1/8" (3.2mm) and 2 threaded barbed fittings (6 to 9mm)	S90	SW90-2-2L-FSC
	Metallic barrel	SWBA-2-2L-FSC
	Space-saving can	SWB-2-2L-FSC
	S55	SW55-2-3L-FSC
Funnel with hinged lid with 1	S60	SW60-2-3L-FSC
activated charcoal cartridge, 2 inlets 1/8" (3,2mm) and 3	S90	SW90-2-3L-FSC
threaded barbed fittings (6 to 9mm)	Metallic barrel	SWBA-2-3L-FSC
	Space-saving can	SWB-2-3L-FSC



*Activated charcoal cartridge sold separately (see page 24)






FUNNEL WITH STOPCOCK

It secures the storage of your residual solvents, equipped with an activated charcoal cartridge for filtration of unobstructed vapors. It is possible to add entries and / or threaded barbed fittings by creating your own reference (see below).

For cans	REFERENCE
S51/55	SW55-C-FS
S60/61	SW60-C-FS







<u>Note</u>: Standard references are not equipped with additional inlets. For another configuration, follow the explanation above and create your own reference!



FUNNEL WITH BALL VALVE WITH HIGH SPILL CAPACITY

Funnel with ball valve allowing to quickly pour used solvents, for plastic can from 5 to 50L.



DESCRIPTION	REFERENCE WITHOUT CARTRIDGE	REFERENCE WITH CARTRIDGE *
Funnel with ball valve for can S55	ET55-B	ET55-C-B
Funnel with ball valve for can S60	ET60-B	ET60-C-B
Funnel with ball valve for can S90	ET90-B	ET90-C-B
Funnel with ball valve for space-saving can S50	ETB-B	ETB-C-B
Funnel with ball valve for metallic barrel	ETBA-B	ETBA-C-B

^{*} Activated charcoal cartridge are to be ordered separately

FUNNEL WITH BALL VALVE WITH SMART WASTE CAP FUNCTION

Funnel with ball valve for 5 to 20L can and metallic barrel, with or without external solvent inlet. It can be equipped with the same functions as a conventional Smart Waste Cap (threaded barbed fittings, activated charcoal cartridge, inlets) to directly receive the effluent from your HPLC chains without passing through the funnel.





Funnels with ball valve with activated charcoal cartridge



For any other configuration, contact us directly!

DESCRIPTION	THREAD	REFERENCE
Funnel with ball valve with an adapter for the activated charcoal cartridge, 2 fittings 1/8" and 1 threaded barbed fitting (6 to 9 mm)	S55	ET55-2-C-1L
	S60	ET60-2-C-1L
	S90	ET90-2-C-1L
	NALGENE® B83	ETB83-2-C-1L

FUNNEL WITH BALL VALVE WITH HIGH SPILL CAPACITY AND CARTRIDGE OF ACTIVATED CHARCOAL FOR METALLIC BARREL

DESCRIPTION	REFERENCE
Funnel with ball valve for metallic or plastic barrel	ЕТВА-В
Funnel with ball valve for metallic or plastic barrel + 1 threaded barbed fitting (6 to 9 mm)	ETBA-1L
Funnel with ball valve for metallic or plastic barrel + 2 threaded barbed fittings (6 to 9 mm)	ETBA-2L
Activated charcoal cartridge of 250 gr for a metallic barrel of 50 L Thread*: 2" for the large bung and 3/4"for the small bung	FC-250
Activated charcoal cartridge of 500 gr for a metallic barrel of 200 L Thread*: 2" for the large bung and 3/4"for the small bung	FC-500

^{*} For other types of thread on your barrel, please contact us directly.





REFERENCE
TAM-PE *
TAM-LUX *

^{*} Add the option reference at the end of the funnel reference.





ELECTRO CONDUCTIVE FUNNEL





DESCRIPTION	THREAD	Reference
Electro conductive funnel with ball valve equipped with a stainless steel sieve	S55	ET55-B-EL
	S60	ET60-B-EL
	S90	ET90-B-EL
	NALGENE® B83	ETB83-B-EL
Electro conductive funnel with ball valve equipped with a stainless steel sieve and an adapter for an activated charcoal cartridge	S55	ET55-C-EL
	S60	ET60-C-EL
	S 90	ET90-C-EL
	NALGENE® B83	ETB83-C-EL

ELECTRO CONDUCTIVE CANS



VOLUME (LITER)	THREAD	REFERENCE
10	S60 / DIN 61	B10-61-EL
20	S60 / DIN 61	B20-61-EL
30	S60 / DIN 61	B30-61-EL
	0007 2111 01	D00 01 EE